



Parts are produced by polypropylene in grades approved for food packaging & human consumption according to New Code for Foodstuff and beverages 1998, Chapter II, Art.26.

- Decoration Methods for the Tub: IML Self Adhesive Label Dry Offset Printing Shrink Sleeve Label
Decoration Methods for the Lid: IML Self Adhesive Label Dry Offset Printing Shrink Sleeve Label
Content's Safety can be Provided By: Tamper Evident Container Heat Sealing Tamper Sleeve

Unit weight of the tub is 29 gr / Unit weight of the lid is 9 gr / Total unit weight is 38 gr
Filling volume is 1080 ml

ARE TUBS PACKED TOGETHER WITH LIDS IN THE SAME CARTON BOX? Yes No

Carton Box Dimensions (for TUBs): 40*65*54 cm Quantity / Box (for TUBs): 375 pieces
Pallet Dimensions (for TUBs): 80*120*230 cm Quantity / Pallet (for TUBs): 7.500 pieces
Carton Box Dimensions (for LIDs): 40*62*34 cm Quantity / Box (for LIDs): 1.125 pieces
Pallet Dimensions (for LIDs): 80*120*253 cm Quantity / Pallet (for LIDs): 31.500 pieces

Loading Quantities (FCL and FTL)

FCL (20'DV) (With Pallets): 64.500 pieces FCL (40'HC) (With Pallets): 145.300 pieces FTL (With Pallets): 199.900 pieces

Notes:

.....
.....
.....